
# MR ALPHABET SOUP: UNDERSTANDING THE PRINCIPLES OF COMMON MAGNETIC RESONANCE ABBREVIATIONS

**Educational Exhibit** 

## OUTLINE

- Introduction to Magnetic Resonance
  - Current Magnetic Resonance Applications
    - Flavors of MR: MRI, MRA, MRS, NMR etc.
    - Clinical Applications
    - Speaking the Same Language
  - Basics for Understanding MR Physics
    - Magnetic Principles
    - Resonance Principles
    - Physics with Hardware
- Image Acquisition
  - Terms at the Terminal
  - Pulse Sequence Design
    - Spin Echo
    - Gradient Echo
    - Common Clinical Pulse Sequences (FLAIR, DCE, FSE, etc.)
- Tips for Understanding New Terms


## THE ABBREVIATED JUNGLE




# MR PHYSICS: Bo

- All MR systems rely on a main magnetic field to produce a signal.
- The main magnetic field is abbreviated a B<sub>0</sub> and referred to as B zero or B "nought"
- There are primarily three types of magnets used in the clinic:
  - Super Conducting Magnets (0.5T and up)
  - Permanent Magnets (up to 0.3T)
  - Resistive and Electromagnets (up to 0.6T)
- The direction of  $B_0$  can vary depending on the types of system.





**Permanent Magnet** 




## TERMS AT THE TERMINAL

- MRI has the ability to acquire a variety of images in relation to contrast and spatial characteristics.
- The parameters you set at the scanner can drastically change the final image characteristics
- Terms you can tweak at the terminal:
  - TR Repetition Time
  - TE Echo Time
  - TI Inversion Time
  - BW Bandwidth
  - **FOV** Field of View
  - NSA or NEX Number of Signal Averages / Number of Excitations
  - N<sub>PE</sub> Number of Phase Encodes
  - N<sub>FE</sub> Number of Frequency Encodes
- Changes you will see in the images:
  - Noise
  - Contrast to Noise
  - Signal to Noise
  - Scan Time
- By understanding principles behind the terms related to image acquisition you can better predict how images will change during protocol optimization

# THE BUILDING BLOCKS OF PSD

- Essentially all pulse sequences can be group into two categories: spin echo (SE) and gradient echo (GE or GRE)
- Both pulse sequences use a radio frequency (RF) pulse to tip the net magnetization into the transverse plane
- SE applies a 180° refocusing pulse to give the image a T2 weighting while GE does not use a refocusing pulse and will have a T2/T2\* weighting.
- From these two basic sequences a multitude of new sequences are born.



# TIPS TO DECODE NEW ABBREVIATIONS

- When learning or using MRI it is common to come across unfamiliar abbreviations.
- Tips for learning a new abbreviation
  - Take a breath and calm down. MRI is confusing and takes time to become familiar with.
  - Determine if the abbreviation is related to physics principles or pulse sequences. This can be made easier by looking for familiar abbreviations contained within the new one. (ex. GRASE → GRASE → SE = Spin Echo)
  - Determine which vendor the term is coming from. (GE,SIEMENS,PHILIPS,ETC.) Similar terms and pulse sequences can change between vendors. (Philips = GRASE, Siemens = TGSE (turbo gradient spin echo)
  - Start from the ground up. Strengthen your understanding of the basics it will help simplifying new complex terms.